

Breast cancer CAD tools

António Polónia, MD, PhD

Ipatimup Diagnostics


@PoloniaAntonio




Disclosure

None


Breast cancer CAD tools



a) normal tissue




b) benign lesion




c) carcinoma *in situ*




d) invasive carcinoma



Breast cancer CAD tools



A



B




Fig. 3. Example of a pixel-wise annotated whole-slide image from the training set.
■ benign; ■ *in situ*; ■ invasive.

Breast cancer CAD tools

Table 2

Summary of the methods submitted for Part A. ^Adetailed description in Section 3.2; ^{AB}detailed description in Section 3.4. Pre-training is performed on ImageNet (Krizhevsky et al., 2012). **Acc.** is the overall prediction success for the four classes; **Approach** lists the main methods to label the images; **Ensemble** (Ens.) indicates if the approach uses a single or multiple models (and their number, when available); **External sets** indicates if the method was trained using datasets other than from Part A; **Context (area ratio)** is the ratio between the original size and the size of the patch used for training the network (prior to rescaling); **Input size (pixels)** is the size of the image to be analyzed by the model; **Color normalization** (color norm.) indicates if any histology-inspired normalization was used. N/A: information not available. ¹Pre-trained on CAMELYON (<https://camelyon17.grand-challenge.org/>). M1: (Bejnordi et al., 2016); M2: (Krishnan and Shah, 2012); M3: (Macenko et al., 2009); M4: (Reinhard et al., 2001).

Team	Acc.	Approach	Pre-trained	Ens.	External sets	Context (area ratio)	Input size (pixels)	Color norm.
(Chennamsetty et al., 2018) (216) ^A	0.87	Resnet-101; Densenet-161	✓	3	✗	1	224 × 224	✗
(Kwok, 2018) (248) ^{AB}	0.87	Inception-Resnet-v2	✓	3	Part B	0.71	299 × 299	✗
(Brancati et al., 2018) (1) ^A	0.86	Resnet-34, 50, 101	✓	3	✗	1	308 × 308 615 × 615	✗
(Marami et al., 2018) (16) ^{AB}	0.84	Inception-v3	✓	4	Part B BreakHis	0.33	512 × 512	M1
(Kohl et al., 2018) (54) ^{AB}	0.83	Densenet-161	✓ ¹	✗	✗	1	205 × 154	✗
(Wang et al., 2018a) (157) ^A	0.83	VGG16	✓	✗	✗	0.765	224 × 224	✗
Steinfeldt et al. (186)	0.81	Xception	✓	✗	✗	0.028-0.751	229 × 229	✗
(Koné and Boulmane, 2018) (19) ^A	0.81	ResNeXt50	✓	✓	BISQUE	1	299 × 299	✗
Nedjar et al. (36)	0.81	Inception-v3, Resnet-50, MobileNet	✓	✓	✗	1	224 × 224	✗
Ravi et al. (412)	0.8	Resnet-152	✓	✗	✗	0.875	224 × 224	M2
(Wang et al., 2018b) (22)	0.79	VGG16	✓	✗	✗	0.255	224 × 224	M3
(Cao et al., 2018) (425)	0.79	PTFAS+GLCM, ResNet-18, ResNeXt, NASNet-A, ResNet-152, VGG16, Random Forest SVM	✓	✓	✗	1	224 × 224 331 × 331	✗
Seo et al. (60)	0.79	ResNet, Inception-V3, Random Forests	✓	✗	✓	1	299 × 299	✗
Sidhom et al. (370)	0.78	ResNet-50	✓	✗	✗	0.018 0.289	224 × 224	M4
(Guo et al., 2018) (242)	0.77	GoogLeNet	✓	2	✗	1 0.083	224 × 224	M4
Ranjan et al. (61)	0.77	AlexNet	✓	2	✗	1	224 × 224	✗
(Mahbod et al., 2018) (73)	0.77	ResNet-50, ResNet-101	✓	2	✗	1	224 × 224	M3
(Ferreira et al., 2018) (18)	0.76	Inception-ResNet-v2	✓	✗	✗	1	224 × 224	✗
(Pimkin et al., 2018) (256)	0.76	ResNet34, Densenet169, Densenet201 XGBoost	✓	12	Part B BreakHis	1	300 × 300	✗
Sarker et al. (358)	0.75	Inception-v4	✓	✗	✗	0.083	299 × 299	✗
(Rakhiin et al., 2018) (98)	0.74	VGG16, ResNet-50, lceptionV3, LightGBM	✓	✓	✗	0.20 0.54	400 × 400 600 × 600	M3
(Iesmantas and Alzbutas, 2018) (164)	0.72	Custom CNN (Capsule Network)	✗	✗	✗	0.029	512 × 512	M4
Xie et al. (253)	0.72	CNN	✗	✗	✗	0.083	512 × 512	✗
(Weiss et al., 2018) (268)	0.72	Xception, Logistic Regression	✓	✗	✗	1	1024 × 768	M3
(Awan et al., 2018) (6)	0.71	ResNet50, SVM	✓	✗	✗	0.33	512 × 512	M4
Liang (62)	0.7	VGG16, VGG19, ResNet50, InceptionV3, Inception-Resnet, k-NN	✓	5	✗	0.083	N/A	✗

Breast cancer CAD tools



Table 3

Summary of the methods submitted for whole-slide image analysis (Part B). ^Bdetailed description in Section 3.3; ^{AB}detailed description in Section 3.4. Pre-training is performed on ImageNet (Krizhevsky et al., 2012) unless stated otherwise. **Score** is the custom metric from Eq. (1); **Approach** lists the main methods to label the images; **Ensemble** (Ens.) indicates if the approach uses a single or multiple models (and their number, when available); **External sets** indicates if the method was trained using datasets other than from Part A; **Context (area ratio)** is the ratio between the average original size and the size of the patch that is used for training the network (prior to rescaling); **Input size (pixels)** is the size of the image to be analyzed by the model; **Color normalization** (Color norm.) indicates if any histology-inspired normalization was used. ¹ trained on ImageNet, ² trained on VOC2012.

Team	Score	Approach	Pretrained	Ens.	External sets	Context (area ratio)	Input size (pixels)	Color norm.
(Kwok, 2018) (248) ^{AB}	0.69	Inception-Resnet-v2	✓	×	Part A	8.5e−4	299 × 299	×
(Marami et al., 2018) (16) ^{AB}	0.55	Inception-v3 + adaptive pooling	✓	4	Part A	9.9e−5	512 × 512	×
Jia et al. (296)	0.52	ResNet-50 + multiscale atrous convolution	✓	×	BreakHis	9.9e−5	512 × 512	×
Li et al. (137)	0.52	VGG16 ¹ , DeepLabV2 ¹ , Resnet50 ²	✓	✓	×	9.9e−5	512 × 512	×
Murata et al. (91)	0.50	U-Net	×	×	×	1.6e−2	256 × 256	×
(Galal and Sanchez-Freire, 2018) (264) ^B	0.50	DenseNet	×	×	×	1.6e−3	2048 × 2048	×
(Vu et al., 2018) (166) ^{AB}	0.49	DenseNet, SENet, ResNext	✓	×	×	1.5e−4	630 × 630	×
(Kohl et al., 2018) (54) ^{AB}	0.42	Densenet-161	✓	×	Part B non-annotated	9.3e−6	157 × 157	×

Breast cancer CAD tools



**Multiclass Classification of Breast Cancer
in Whole-Slide Images**

Scotty Kwok^(✉)
Seek AI Limited, Hong Kong, China





Fig. 1. Overview of the framework

@PoloniaAntonio

Breast cancer CAD tools



Breast cancer CAD tools



Breast cancer CAD tools




Table 3

Diagnostic Accuracy in Test A and B

	Phase	Test A	Test B
Algorithm A	1	0.87	NA
Algorithm B	1	NA	0.49
Pathologists (average)	1	0.83	0.83
	2	0.87	0.82
	3	0.90	NA
Residents (average)	1	0.77	0.89
	2	0.82	0.88
	3	0.87	NA
All observers (average)	1	0.80 ^a	0.86
	2	0.85 ^{a,b}	0.85
	3	0.88 ^b	NA

NA, not applicable.

^aWilcoxon, $P < .001$

^bWilcoxon, $P = .001$.

Breast cancer CAD tools




Table 3

Diagnostic Accuracy in Test A and B

	Phase	Test A	Test B
Algorithm A	1	0.87	NA
Algorithm B	1	NA	0.49
Pathologists (average)	1	0.83	0.83
	2	0.87	0.82
	3	0.90	NA
Residents (average)	1	0.77	0.89
	2	0.82	0.88
	3	0.87	NA
All observers (average)	1	0.80 ^a	0.86
	2	0.85 ^{a,b}	0.85
	3	0.88 ^b	NA

NA, not applicable.

^aWilcoxon, $P < .001$

^bWilcoxon, $P = .001$.

Breast cancer CAD tools




Table 3

Diagnostic Accuracy in Test A and B

	Phase	Test A	Test B
Algorithm A	1	0.87	NA
Algorithm B	1	NA	0.49
Pathologists (average)	1	0.83	0.83
	2	0.87	0.82
	3	0.90	NA
Residents (average)	1	0.77	0.89
	2	0.82	0.88
	3	0.87	NA
All observers (average)	1	0.80 ^a	0.86
	2	0.85 ^{a,b}	0.85
	3	0.88 ^b	NA

NA, not applicable.

^aWilcoxon, $P < .001$

^bWilcoxon, $P = .001$.

Breast cancer CAD tools




Table 3

Diagnostic Accuracy in Test A and B

	Phase	Test A	Test B
Algorithm A	1	0.87	NA
Algorithm B	1	NA	0.49
Pathologists (average)	1	0.83	0.83
	2	0.87	0.82
	3	0.90	NA
Residents (average)	1	0.77	0.89
	2	0.82	0.88
	3	0.87	NA
All observers (average)	1	0.80 ^a	0.86
	2	0.85 ^{a,b}	0.85
	3	0.88 ^b	NA

NA, not applicable.

^aWilcoxon, $P < .001$

^bWilcoxon, $P = .001$.

Breast cancer CAD tools




Table 3

Diagnostic Accuracy in Test A and B


	Phase	Test A	Test B
Algorithm A	1	0.87	NA
Algorithm B	1	NA	0.49
Pathologists (average)	1	0.83	0.83
	2	0.87	0.82
	3	0.90	NA
Residents (average)	1	0.77	0.89
	2	0.82	0.88
	3	0.87	NA
All observers (average)	1	0.80 ^a	0.86
	2	0.85 ^{a,b}	0.85
	3	0.88 ^b	NA

NA, not applicable.


^aWilcoxon, $P < .001$

^bWilcoxon, $P = .001$.


Breast cancer CAD tools



Breast cancer CAD tools



Breast cancer CAD tools



A1 6/7 correct
A3 5/7 correct


A1 5/7 correct
A3 no change

A1 3/7 correct
A3 2/7 correct


B1 6/6 correct
B2 no change

Output of the AI: invasive carcinoma (incorrect)


Breast cancer CAD tools



A1 3/7 correct
A3 5/7 correct




A1 4/7 correct
A3 6/7 correct



A1 3/7 correct
A3 5/7 correct

Output of the AI: carcinoma in situ (correct)

Breast cancer CAD tools




HE

CK5/6


ER

Output of the AI: carcinoma in situ (correct)

Breast cancer CAD tools




A1 1/7 correct
A3 1/7 correct



B1 5/6 correct
B2 4/6 correct


Output of the AI: normal (incorrect)

Breast cancer CAD tools




A1 1/7 correct
A3 1/7 correct

Output of the AI: normal (incorrect)




E-cadherin

Breast cancer CAD tools



A1 0/7 correct
A3 1/7 correct




B1 2/6 correct
B2 3/6 correct


Output of the AI: benign (correct)

Breast cancer CAD tools


HE




CK5/6



p63



ER



Output of the AI: benign (correct)

Breast cancer CAD tools

Most clinically relevant incorrect classifications

Benign vs Carcinoma in situ

Test A

Phase 1 - 7,1%

Phase 3 - 3,4%



Test B


Phase 1 - 6,3%

Phase 2 - 6,5%




Breast cancer CAD tools


Concordance with AI
 $P=0.018$ (Wilcoxon)



Breast cancer CAD tools




Breast cancer CAD tools




Breast cancer CAD tools

Probability of being correct if discordant



Phase 1
P=0.091 (Wilcoxon)




Phase 3
P=0.310 (Wilcoxon)

@PoloniaAntonio


Breast cancer CAD tools

Probability of being correct if discordant (Experts)
 $P=0.018$ (Wilcoxon)




Breast cancer CAD tools

Probability of being correct (Experts)
P=0.018 (Wilcoxon)



Breast cancer CAD tools

Will machines replace humans?



Breast cancer CAD tools

Machines need the Holy Grail!!!



The Holy Grail is the Ground Truth!



Breast cancer CAD tools

It's TIME to choose – red pill or blue pill!



Breast cancer CAD tools

